# Package 'RLumSTARR'

May 10, 2022
Type Package
Title Spatially Resolved Radiofluorescence Analysis (EXPERIMENTAL PACKAGE)
Version 0.1.0.9000-132
Author Sebastian Kreutzer [aut, cre] ( <a href="https://orcid.org/0000-0002-0734-2199">https://orcid.org/0000-0002-0734-2199</a> ),  Dirk Mittelstrass [aut] ( <a href="https://orcid.org/0000-0002-9567-8791">https://orcid.org/0000-0002-9567-8791</a> )
Maintainer Sebastian Kreutzer < sebastian.kreutzer@uni-heidelberg.de>
Description Analysing spatially resolved radiofluorescence data using ImageJ in combination with Bayesian number-crunching.  A collection of functions to support early work on the subject.
License GPL-3
<b>Depends</b> R (>= 4.0), utils,
Imports cli (>= 3.2.0),
Suggests rmarkdown (>= 2.11), rjags (>= 4-12), testthat (>= 3.1.2), knitr (>= 1.37)
Encoding UTF-8
Language en-GB
LazyData true
VignetteBuilder knitr
RoxygenNote 7.1.2
R topics documented:
RLumSTARR-package

2 create\_RFCurveArray

	get_MCMCParame melt_RLumSTARR																				7
	run_ImageJ run_TRUELightExt																				
Index	_ 0																				12
RLumS	TARR-package	SpaTiA	Ally	Re	sol	ved	Rι	ıdi	ofl	luo	re:	sce	enc	ce							

## **Description**

A collection of functions to analyse spatially resolved radiofluorescence data

#### **Details**

## **Funding**

• Sebastian Kreutzer received funding from the European Union's Horizon 2020 research and innovation programme under the Marie Skłodowska-Curie grant agreement No 844457 (CREDit).

## Author(s)

Sebastian Kreutzer, Institute of Geography, Heidelberg University (Germany) Dirk Mittelstrass (Germany), RLum Developer Team

## References

##TODO

create\_RFCurveArray Create Multidimensional Curve Arrays from RF File input

## **Description**

Helper function to create a multidimensional curve array based on RF-file input imported using the function Luminescence::read\_RF2R() to prepare the Bayesian modelling process

## Usage

```
create_RFCurveArray(files)
```

## **Arguments**

files

list (**required**): list of .rf files to be imported. Alternatively you can pass an object created by Luminescence::extract\_ROI. In this case only the first list element RF\_nat is filled.

## Value

Returns a list of class RLumSTARR\_RFCurveArrary with two arrays for the RF\_nat and the RF\_reg curve

extract\_TRUELight 3

#### **Function version**

0.1.0

#### How to cite

Kreutzer, S., 2022. create\_RFCurveArray(): Create Multidimensional Curve Arrays from RF File input. Function version 0.1.0. In: Kreutzer, S., Mittelstrass, D., 2022. RLumSTARR: Spatially Resolved Radiofluorescence Analysis (EXPERIMENTAL PACKAGE). R package version 0.1.0.9000-132.

#### Author(s)

Sebastian Kreutzer, Institute of Geography, Heidelberg University (Germany) , RLum Developer Team

#### See Also

Luminescence::read RF2R, Luminescence::extract ROI

## **Examples**

```
## list files using package external data
files <- list.files(system.file("extdata", "", package="RLumSTARR"), full.names=TRUE)
## create curve array
create_RFCurveArray(files = files)</pre>
```

extract\_TRUELight

Extract True Light from the Camera Measurements using a Bayesian Approach

## Description

A Bayesian modelling approach to extract the true light using the expanding region-of-interest (ROI) approach proposed by Cunningham and Clark-Balzan (2017). The function will return the results for a **single** curve from a **single** ROI

#### Usage

```
extract_TRUELight(
  data,
  element = c("RF_nat", "RF_reg"),
  ROI = 2,
  stepping = 1,
  method_control = list(),
  verbose = TRUE
)
```

4 extract\_TRUELight

#### **Arguments**

data array (required): object created by create\_RFCurveArray

element character (with default): element from the input to be analysed, 'supported are

only RF\_nat or RF\_reg

ROI numeric (optional): ROI to be analysed, if nothing is given all ROIs are anal-

ysed, however, the first ROIS is discarded!

stepping numeric (with default): stepping parameter that allows you to model only every

xth (the value in stepping). This option can be extremely useful to play with data because it dramatically improves the modelling speed because less data are

considered.

method\_control list (optional): parameter to be passed to rjags. Supported are n.chain, n.iter,

thin, variable.names, model, see details for more.

verbose logical (with default): enable/disable terminal feedback

#### **Details**

#### Method control

Supported options to be passed via the parameter method\_control, most of them are used internally for the calls to runjags::run.jags and rjags::coda.samples.

#### **PARAMETER**

n.chain

thin

burnin

sample

adapt

summarise

method

jags.refresh\tab [numeric] \tab0.1\tab refresh rate of update of the iteration process, select larger

Note: The argument model allows to heavily modified the underlying model. To avoid crashes the paramters passed by variable.names will always be cross-checked against parameteres present in the model. Unknown parameters will be skipped!

#### Value

Returns a list of class RLumSTARR\_TRUELight with an the following elements:

- ....\$RF\_curve: Luminescence::RLum.Data.Curve object (the RF curve with the true light)
- ...\$rjags\_output: rjags::coda.samples output for further processing. Note: Regardless the observed variable, the parameter alpha will always be used to create the curve
- $\dots$ \$model: the model used to run the Bayesian process, use writeLines to have nicely formatted terminal output

#### **Function version**

0.1.1

get\_MCMCParameters 5

#### How to cite

Kreutzer, S., 2022. extract\_TRUELight(): Extract True Light from the Camera Measurements using a Bayesian Approach. Function version 0.1.1. In: Kreutzer, S., Mittelstrass, D., 2022. RLum-STARR: Spatially Resolved Radiofluorescence Analysis (EXPERIMENTAL PACKAGE). R package version 0.1.0.9000-132.

#### Author(s)

Sebastian Kreutzer, Institute of Geography, Heidelberg University (Germany) , RLum Developer Team

#### References

Cunningham, A.C., Clark-Balzan, L., 2017. Overcoming crosstalk in luminescence images of mineral grains. Radiation Measurements 106, 498–505. doi:10.1016/j.radmeas.2017.06.004

## See Also

```
create_RFCurveArray, get_MCMCParameters
```

## **Examples**

```
## list files using package external data
files <- list.files(system.file("extdata", "", package="RLumSTARR"), full.names=TRUE)

## create curve array
dat <- create_RFCurveArray(files = files)
output <- extract_TRUELight(
    data = dat,
    ROI = c(4),
    stepping = 10,
    verbose = FALSE,
    method_control = list(
        n.chain = 2,
        sample = 100,
        thin = 20))</pre>
```

get\_MCMCParameters

Extracts a parameter from an MCMC list

## Description

Short cut to extract a parameter from an MCMC list or an object of class RLumSTARR\_TRUELight created by the functions extract\_TRUELight or run\_TRUELightExtraction. In the latter case, the row names of the extracted matrices have the same dimension as the RF curve. If more processing is wanted, 'coda' package should be used.

#### Usage

```
get_MCMCParameters(mcmc, parameter = NULL, prob = 0.95, unlist = TRUE)
```

get\_MCMCParameters

#### **Arguments**

mcmc coda::mcmc or coda::mcmc.list (required): input object, if created by extract\_TRUELight

the correct object is extracted automatically

parameter character (optional): name of the parameter to be extracted. If NULL (the default)

all found parameters are extracted and the result is a list with matrices of those parameter. If parameter is a vector, the function will try to extract the names

parameters.

prob numeric (with default): probability for the HPD calculation (cf. coda::HPDinterval)

unlist logical (with default): if TRUE the output is a matrix of the means of the lower

and upper intervals of the parameter. If the parameter was estimated based on multiple chains, this chains are also subject to an average calculation. If FALSE

the output is a list as returned by coda::HPDinterval

#### Value

Returns a matrix with the parameter value or a named list with such matrices if parameters has a length > 1.

#### **Function version**

0.1.1

## How to cite

Kreutzer, S., 2022. get\_MCMCParameters(): Extracts a parameter from an MCMC list. Function version 0.1.1. In: Kreutzer, S., Mittelstrass, D., 2022. RLumSTARR: Spatially Resolved Radiofluorescence Analysis (EXPERIMENTAL PACKAGE). R package version 0.1.0.9000-132.

## Author(s)

Sebastian Kreutzer, Institute of Geography, Heidelberg University (Germany) , RLum Developer Team

#### See Also

```
coda::HPDinterval, extract_TRUELight
```

## **Examples**

```
## load example files
files <- list.files(system.file("extdata", "", package="RLumSTARR"), full.names=TRUE)

##prepare data and run model
dat <- create_RFCurveArray(files = files)
output <-
extract_TRUELight(
   data = dat,
   ROI = c(4),
   stepping = 60,
   verbose = FALSE,
   method_control = list(
   n.chain = 1,
   sample = 50,
   thin = 20))</pre>
```

melt\_RLumSTARR 7

```
##extract parameters
get_MCMCParameters(output)
```

melt\_RLumSTARR

Melt output from RLumSTARR into simple data frames

## Description

The function provides a convenient way to convert the objects created by RLumSTARR into simple data.frames that can be processed conveniently by other functions for example the package ggplot2.

## Usage

```
melt_RLumSTARR(x, ...)
```

## **Arguments**

x (**required**): input object of class RLumSTARR\_RFCurveArray or RLumSTARR\_TRUELight
... further parameters, currently not used

## Value

data.frame

#### **Function version**

0.1.0

## How to cite

Kreutzer, S., 2022. melt\_RLumSTARR(): Melt output from RLumSTARR into simple data frames. Function version 0.1.0. In: Kreutzer, S., Mittelstrass, D., 2022. RLumSTARR: Spatially Resolved Radiofluorescence Analysis (EXPERIMENTAL PACKAGE). R package version 0.1.0.9000-132.

#### Author(s)

Sebastian Kreutzer, Geography & Earth Sciences, Aberystwyth University (United Kingdom) , RLum Developer Team

8 run\_ImageJ

run\_ImageJ

Run ImageJ SR-RF macro

#### **Description**

The script runs the SR-RF ImageJ macro in batch mode out of R

#### Usage

```
run_ImageJ(
 path,
 RF_nat = "default",
 RF_reg = "default",
 bg_rm = "take_from_RF_reg",
  image_group_size = 5,
  image_alignment = TRUE;
  first_slices_rm = FALSE,
 noise_tolorance = 10,
 ROI_size = 10,
  center_x = 0.5,
 center_y = 0.5,
 diameter = 0.9,
  use_predefined_ROIs = FALSE,
  channel_time = 5,
  save_workflow_images = FALSE,
  save_additional_results = FALSE,
  save_signal_decay_videos = FALSE,
 offset_time = 0,
  .ImageJ = "/Applications/Fiji.app/Contents/MacOS/ImageJ-macosx"
)
```

## **Arguments**

```
path
                  character (required): path to files to be analysed
                  character (default): name of the RF_nat file
RF_nat
                  character (default): name of the RF_reg file
RF_reg
                  character (with default): background subtraction options. Allowed are none (no
bg_rm
                  background subtraction), take_from_RF_reg (takes the last 100 channels from
                  the RF_reg signal: dangerous) or <your file name> (this does not work in
                  batch mode)
image_group_size
                  numeric (with default): grouping value for running median to remove outliers
image_alignment
                  logical (with default): enable/disable image alignment
first_slices_rm
                  logical (with default): remove first slice of each curve set
noise_tolorance
                  numeric (with default): noise tolerance parameter
                  numeric (with default): ROI size in pixel
ROI_size
```

run\_ImageJ 9

center\_x numeric (with default): aliquot ROI centre x-coordinate numeric (with default): aliquot ROI centre y-coordinate center\_y diameter numeric (with default): relative diameter aliquot ROI use\_predefined\_ROIs logical (with default): use pre-defined ROIs imported from a file ROIs.zip found in the same folder as the files numeric (with default): channel time, this parameter was set the moment the channel\_time sequence was written save\_workflow\_images logical (with default): enable/disable writing of additional workflow images save\_additional\_results logical (with default): enable/disable writing of additional workflow images save\_signal\_decay\_videos logical (with default): enable/disable writing of additional workflow videos  $offset\_time$ numeric (with default): offset time for the time axis

numeric (with default): Path to ImageJ (the macro is shipped with the package)

#### Value

This functions returns the path of the analysed data

#### **Function version**

.ImageJ

0.1.0

## How to cite

Kreutzer, S., 2022. run\_ImageJ(): Run ImageJ SR-RF macro. Function version 0.1.0. In: Kreutzer, S., Mittelstrass, D., 2022. RLumSTARR: Spatially Resolved Radiofluorescence Analysis (EXPERIMENTAL PACKAGE). R package version 0.1.0.9000-132.

## Author(s)

Sebastian Kreutzer, Institute of Geography, Heidelberg University (Germany) , RLum Developer Team

## **Examples**

##TODO

```
run_TRUELightExtraction
```

Run Automated TRUE-Light Extraction over multiple datasets

#### **Description**

Runs extract\_TRUELight automatically in a parallel processing mode using parallel::mclapply. In essence, the function calls extract\_TRUELight and extracts the RF curves from the results and constructs an Luminescence::RLum.Analysis for further analyses.

## Usage

```
run_TRUELightExtraction(
  data,
  ROI,
  stepping = 1,
  mc.cores = max(c(1, parallel::detectCores() - 2)),
  method_control = list(),
  include_jags_output = FALSE,
  verbose = TRUE
)
```

## Arguments

data array (**required**): object created by create\_RFCurveArray

ROI numeric (optional): ROIs to be analysed, if nothing is given all ROIs are anal-

ysed, however, the first ROIS is discarded!

stepping numeric (with default): the stepping parameter from and to be passed to ex-

tract\_TRUELight

mc.cores numeric (with default): number of cores used for the processing, passed to par-

allel::mclapply

method\_control list (optional): parameters to be passed to extract TRUELight

include\_jags\_output

logical (with default): allows to include the output from JAGS as info object in the Luminescence::RLum.Analysis output objects for further diagnostics. Setting this option to TRUE is not recommended for large datasets, since it will tremendously inflate the size of output and consume a lot of memory. If in-depth diagnostics are required, the function extract\_TRUELight is recommended in-

stead.

verbose logical (with default): enable/disable verbose mode. The output of the MCMC

sampling using rjags is always silent.

## Value

The output is a list containing Luminescence::RLum.Analysis objects with two Luminescence::RLum.Data.Curve objects for RF\_nat and RF\_reg respectively.

## **Function version**

0.1.0

#### How to cite

Kreutzer, S., 2022. run\_TRUELightExtraction(): Run Automated TRUE-Light Extraction over multiple datasets. Function version 0.1.0. In: Kreutzer, S., Mittelstrass, D., 2022. RLumSTARR: Spatially Resolved Radiofluorescence Analysis (EXPERIMENTAL PACKAGE). R package version 0.1.0.9000-132.

## Author(s)

Sebastian Kreutzer, Institute of Geography, Heidelberg University (Germany) , RLum Developer Team

## **Examples**

```
## list files using package external data
files <- list.files(system.file("extdata", "", package="RLumSTARR"), full.names=TRUE)
## create curve array
dat <- create_RFCurveArray(files = files)
output <- run_TRUELightExtraction(
data = dat,
stepping = 15,
mc.cores = 1,
ROI = 5,
verbose = TRUE,
method_control = list(
n.chain = 1,
sample = 0,
thin = 20))</pre>
```

## **Index**

```
* datagen
    create_RFCurveArray, 2
    extract_TRUELight, 3
* package
    RLumSTARR-package, 2
array, 4, 10
character, 4, 6, 8
coda::HPDinterval, 6
coda::mcmc, 6
coda::mcmc.list,6
create_RFCurveArray, 2, 4, 5, 10
data.frame, 7
extract_TRUELight, 3, 5, 6, 10
get_MCMCParameters, 5, 5
list, 2, 4, 6, 10
logical, 4, 6, 8-10
Luminescence::extract_ROI, 2, 3
Luminescence::read_RF2R, 3
Luminescence::RLum.Analysis, 10
Luminescence::RLum.Data.Curve, 4, 10
matrix, 6
melt_RLumSTARR, 7
numeric, 4, 6, 8-10
parallel::mclapply, 10
rjags::coda.samples,4
RLumSTARR (RLumSTARR-package), 2
RLumSTARR-package, 2
run_ImageJ, 8
run_TRUELightExtraction, 5, 10
runjags::run.jags,4
writeLines, 4
```